

Quick user Barcode login

Wednesday, 27 April 2011 19:27 - Last Updated Wednesday, 27 April 2011 20:31



The DRX500Sx scales are equipped with user login and identification features which allow the server to monitor user performed pill counting transactions. Login to the scale is fast and simple. Each scale is supplied with user access cards which are assigned to individual users. To login a user simply scans a barcode on the access card. In large settings such as hospitals or clinics, users can use barcodes printed on employee identification cards. If the scale's login type is set to password a user can also login with a user name and password.

< [Previous](#) | [Next](#) >