| DRX-SSX SYSTEM ** | | | | | | 365///// | Congel via abre gain |
|--|--|---|----|---------|------------|---|----------------------|
| The management of the second s | Minut sting | | | | 22/1/1/1/1 | | |
| States Same Salaria | Owner tasket | Awres | | - | - | 41111111 | |
| | Dente tate | Texts, sta | | - | - | | |
| | DHOP mode: IP Address: And manifest Transport protocoli | materia | 10 | - | - | -//// | Tast. 1 Robert 7 |
| | | 10,00.00 | - | - | - | | |
| | | 1001 | | - | - | m / / / / / / / / | |
| | | No. | 10 | - | - | 2000000 | |
| | Connection Seasouth | | - | - | | 00000000 | |
| | Routing model | alard . | 10 | - | | 00000000 | |
| | Computer made | manufacture | | - | | 10///////// | |
| | Remote hand #1 | 10,10,10,1 | - | - | - | 200000000 | |
| | Remote hand part | 128 | | - | - | 5000000 | |
| | number: Gebrung: Reformit: | | | | | 20000000 | |
| | | 10,08108. | | - | - | 000000 | |
| | | 20.20.2012 | | - | - | | |
| | Percent | | | - | - | 900000000 | |
| | | | | | | 0000000 | hander a to |
| a Parata | Second Second | | | | _ | 00000000 | Charles and the |
| Inference in a second sec | baut de | 1410 | | - | - | 0////////////////////////////////////// | /www.hitesteeter.aut |
| Man, Au Digits | Parity: | | | 0000000 | 00000000 | ///www.selatum///* | |
| and the second se | Bits pair bytes | Bella . | | - | - | Langeline al | |
| | | | | | | 27777777 | |
| | Tre gradation until g | han an a | | | | 10//////// | |
| | Pasimum pacient arright | 28 | | - | - | | |
| | Start on any characteric | materia | | - | - | V//////// | |
| | the start characters | dialog . | 12 | - | - | 11111111 | |
| | the stop characters | distant | | | - | 000000000 | |
| | | | | | | 222222222 | |

The software is viewed and controlled through the use of a web browser such as Internet Explorer or Mozilla Firefox. Each server offers administrative flexibility, and can be configured to fit particular pharmacies needs and requirements. For example, some servers may be configured in a very strict manner in order to track and control all users who perform transactions, whereas other servers may be configured more freely where users are not required to log in, and the server is mainly used for sharing drug data and information. The server provides all essential network configuration parameters that allow the server to function securely in even the most complex network environment.



The TORBAL Rx server is designed to handle a network of up to 30 DRX-500s scales. Up to 7 scales can be connected directly to server and more scales can be connected by adding on a hub or a switch or by connecting scales through an existing local area network, as well as, a remote over the internet. For large chains requirements of 30 or more DRX-500s workstations Fulcrum Inc. offers a Windows OS setup version of the TORBAL Rx server that can be installed on a dedication Windows Server PC.

< <u>Previous</u> | <u>Next</u> >